

Data sheet

Laboratory Shakers

LS 700



The photo above is for reference only, may show additional options not included in standard equipment. The real appearance, particularly color and structure of the material may differ from the ones presented in the photo.


Standard equipment:

- microprocessor control of rotation and time
- LCD digital display
- variable speed control: 30 ... 300rpm
- shaking mode: from 1min to 99h or continuous operation

TECHNICAL DATA	
movement	orbital
controller	microprocessor
display	LCD
speed range [rpm]	30...300
accuracy [rpm]	10
amplitude [mm]	5 or 10
max load capacity[kg]	10
shaking mode	1min. ... 99h or continuous operation
dimensions without/with shaking table [mm]	
width	700
height	120/220
depth	420
fits to cooled incubator	ILW 400
nominal power [W]	60
weight with shaking table [kg]	25
ambient temperature [°C]	+10...+40
humidity [%]	up to 70
voltage	230 V 50 Hz
warranty	24 months
manufacturrer	POL-EKO-APARATURA

ACCESSORIES

universal platform




platform for Petri plates shaking



platform for separatory funnels



platform for fixing flasks handles



anti-skid mat

